## [M-Stars D] (For research use only )

Product
Description

M-Stars D (Keratinocyte Growth Medium) is a medium developed for culture of epithelial cells such as cornea, conjunctiva, skin and oral mucosa. M-Stars D is suitable for primary cell culture of corneal epithelial cells and oral mucosal epithelial cells. This medium contains growth factor, and therefore it is possible to use of cell culture immediately by a suitable addition of FBS and antibiotic as the need arises.

**Product** 

M-Stars D 500 mL

Contents

Product No.

Storage

2 to 8°C in the dark

Shelf-life

This medium has a shelf-life of 10 weeks if stored in the dark at 2 to 8°C.

Instructions

for Use

- 1. M-Stars D incubated at 37°C after addition of FBS and antibiotic as the need arises.
- 2. Cells are gently suspended in M-Stars D.
- 3. Cultures are incubated in atmosphere of 5% CO<sub>2</sub> in air at 37°C.

Cell Culture
Conditions

- 1. In case of use amniotic membranes as matrix for cultures, this medium demonstrates superior cell growth of corneal epithelial cells and oral mucosal epithelial cells<sup>1)</sup>.
- 2. Cell growth activates by use of feeder cells.
- 3. Since cell growth has individual difference, it is necessary to examine concentration of FBS.

Caution

- 1. For research use only.
- 2. Not intended for human or animal diagnostic or therapeutic uses.
- 3. In case of low cell density, cell growth may take a lot of time to proliferation.
- 4. Since cell growth has individual difference, in some cases cells may not proliferate.

References

- 1) Nakamura T, Ang LPK, Rigby H, Sekiyama E, Inatomi T, Sotozono C, Fullwood NJ, Kinoshita S. The Use of Autologous Serum in the Development of Corneal and Oral Epithelial Equivalents in Patients with Stevens-Johnson Syndrome. *Invest Ophthalmol Vis Sci.* 2006; 47: 909–916
- For further information on this products, contact at the following.



Kobe International Business Center Bldg.

5-5-2 Minatojima-minamimachi,Chuo-ku,

Kobe-city, Hyogo 650-0047, Japan

TEL. +81-78-304-5870 FAX. + 81-78-304-5871

E-mail: info@arblast.jp